

## **Thiol Disulfide exchange protocol**

### **Extracting yeast mitochondria from the culture (Meisinger, Pfanner, & Truscott, 2006)**

1. Check that samples reach an OD<sub>600nm</sub> of ~0.8
2. Centrifuge 50mL of yeast culture at 3000g for 5 mins a
3. Discard the supernatant
4. Resuspend and wash pellet with ~300mL of pure water
5. Centrifuge at 3000g for 5 mins and discard the supernatant
6. Resuspend pellet with 5mL pre-warmed DTT buffer 100 mM Tris-H<sub>2</sub>SO<sub>4</sub>, pH 9.4, 10 mM dithiothreitol (DTT), prewarmed to 30°C. Add DTT prior to use. A 1 M Tris-H<sub>2</sub>SO<sub>4</sub> stock can be stored at room temperature.
7. Shake yeast and DTT buffer mixture slowly at 30°C for 20 mins
8. Centrifuge at 3000g for 5 mins and discard supernatant
9. Resuspend pellet in 25mL zymolyase buffer
10. Centrifuge at 3000g for 5 mins and discard supernatant
11. Resuspend pellet in 25mL zymolyase buffer
12. Shake solution slowly at 30°C for 30-45 mins
13. Centrifuge at 3000g for 5 mins and discard supernatant
14. Resuspend pellet in 25mL zymolyase buffer
15. Centrifuge at 3000g for 5 mins and discard supernatant
16. Resuspend pellet in 30mL ice-cold homogenisation buffer (0.6 M sorbitol, 10 mM Tris-HCl, pH 7.4, 1 mM ethylenediaminetetraacetic acid (EDTA), 1 mM phenylmethylsulfonyl fluoride (PMSF), 0.2% (w/v) bovine serum albumin (BSA; essentially fatty acid-free, Sigma-Aldrich, Taufkirchen, Germany). Add PMSF from a freshly prepared 100 mM stock in ethanol just prior to use.)
17. Using a Dounce homogeniser, homogenise the spheroplasts with ~15 strokes at 4°C
18. Use 60mL of homogenisation buffer to wash out and collect as much of the sample as possible into a falcon tube
19. Centrifuge homogenate at 1500g for 5 mins at 4°C to pellet cell debris and nuclei. Collect the supernatant
20. Centrifuge the supernatant at 4000g for 5 mins at 4°C. Collect the supernatant and discard the pellet. Collect 50uL aliquots from each sample for 1D gel and quantitation.
21. Centrifuge the supernatant at 12 000g for 15 mins at 4°C. Keep the pellet

22. Resuspend the mitochondrial pellet in SEM (50 mM sucrose, 1 mM EDTA, 10 mM MOPS-KOH, pH 7.2)

23. To store the sample, use liquid nitrogen to flash freeze and store at -80oC